

Lipids

Mission Ultra Cholesterol

SPECIFICATIONS & FEATURES

- 4 μ L Capillary whole blood
- Total CHOL results in as little as 15 seconds
- Biosensor technology for direct sampling application
- Memory for 200 results
- Operated in battery or standard AC adaptor mode



CE 0123



Mission Cholesterol Meter

SPECIFICATIONS & FEATURES

- 10 μ L (individual test), 35 μ L (3-in-1 lipid panel test)
- Results in \leq 2 minutes
- Capillary or venous whole blood, serum or plasma
- Accuracy comparable to laboratory analyzers
- Reflectance Photometry



CE 0123



MISSION CHOLESTEROL 3 IN 1 LIPID TEST 3 IN 1

- Tests for the complete lipid panel
 - Total Cholesterol (CHOL)
 - High-Density Lipoprotein (HDL)
 - Triglycerides (TRIG)
 - Calculated Low-Density Lipoprotein, CHOL/HDL ratio
 - Cardiac Risk Assessment (CHD)
- 18-months extended shelf life
- Room temperature storage
- Individually pouched test devices preserve stability



MISSION CHOLESTEROL METER FEATURES

- Large, easy-to-read display
- Fast, reliable lipid test results
- Compact, ergonomic design
- Memory, stores up to 200 results