



FM 150 Vein trainer

german engineering



Air pressure massage with inflating and deflating air cushion



Continuously adjustable massage intensity



Improves the venous blood circulation



- Revitalising pressure massage
- Promotes blood circulation and relieves heavy and tired legs
- Looses tension
- Air pressure massage with inflating and deflating air cushion
- Continuously adjustable massage intensity
- Alternate inflation of the leg cuffs
- Incl. Safety STOP button
- With deflation function - for emptying the air cushions quickly
- Incl. timer function 10/20/30 min.
- Easy operation with manual control switch
- Versatile application - whether you're at home or on the go
- Individual adjustment of the leg cuff using hook-and-loop fasteners
- Storage pouch included
- Medical device
- Battery or mains adaptor driven
- Mains operation:
- Input: 100-240 V~; 50/60 Hz; 200 mA
- Output: 6 V $\overline{\text{---}}$ 500 mA; 3 W
- Incl. 4 x 1.5V battery AA
- Product measurements (LxWxH): Leg cuff: 61 x 67 x 0.5 cm
- Product weight: approx. 1.4 kg
- 3 year guarantee
- Sales unit: 2 / Shipping carton: -
- EAN no.: 4211125 64325 6
- Item no.: 643.25

643.25_2020-01-14_03_DS_EN_BEU Subject to errors and changes

PhysioLine