

DIGISKOP[®] | VIDEO OTOSCOPE VCLICK30

FEATURES

- High resolution image with 1024x720 pixels.
- Easy to carry hand-held device for moveable diagnosis treatments.
- Wide angle lens provide wide field of view.
- Photo and video recording with screen rotation for easy observation and review.
- Excellent ability to store data with memory of 32GB.
- Ideal device for cerumen management.
- Light weight, compact and ergonomic design.
- Efficient and cost-effective otoscopic examinations.
- Special disposable devices to avoid infection and quickly inserting and removing of scope tube.
- **Outer ear canal:** tympanic membrane, and detection of hearing aids devices.
- USB type-C of the charging base.
- USB port to connect to the PC, MAC or Laptop.
- Charging indicator light.
- Docking charge for easy access.
- Rechargeable lithium battery with capacity of 1350mAh.



REFERENCE

Digital Otoscope Host W/ Aluminum Handgrip, Built-in LCD, No wiring, Built-in battery with Disposable Ear Tips

60.550.110 

